

Specifications for MagentIQ THERAPY AMPF

BASIC CHARACTERISTICS

The device model:	MagentIQ THERAPY AMPF
Number of programs	11+4
Duration of therapy	15 minutes
Frequency magnetic field	1,2 - 9999 Hz
Strength of the magnetic field	0.012 - 50 Gauss
The voltage of the input power	110 - 240 VAC
The frequency of the input power	47 - 63 Hz
Power input	42 VA
Device dimensions (approx.)	8.9x2.8x7.8 inches
GVWR (approx.)	1.98 lbs
OPERATING CONDITIONS AT WORK	
Temperature	5°F to 104°F
Relative humidity	30% to 75%
Atmospheric pressure	20.6 inHg to 31.3 inHg

www.hsttsolutions.com



Available soon in Canada, Qatar, Pakistan, India, Australia and Japan.

magentiQ
THERAPY

Good
vibrations.

ALTERNATE MAGNETIC
PULSATING FIELD

HSTT
SOLUTIONS

SMART CLINIC

**Health Screening, Telemedicine
and Advanced Therapies LLC**

8588 Katy Freeway | Suite 230
Houston, Texas 77024

Office +1.713.234.6861

Fax +1.781.827.0849

www.hsttsolutions.com



HSTT SOLUTIONS

SMART CLINIC

What is MagnetIQ THERAPY AMPF?

MagnetIQ THERAPY Alternate Magnetic Pulsating Field (AMPF) is an innovative device that provides pulsating bio-magnetic resonance therapy. It was developed based on research which respects the physiological processes in the body.

MagnetIQ THERAPY AMPF excites the movement of electrically charged particles (ions) in the cell using a pulsed magnetic field and a specific frequency, so that there is several times higher permeability of the cell membrane. This improves the input of nutrients and oxygen to cells, and the output of toxins and heavy metals.

MagnetIQ THERAPY AMPF has a choice of 11 therapies and 4 programs to stimulate brain waves (alpha, beta, delta and theta). Duration of one session is 15 minutes.

MagnetIQ THERAPY AMPF provides a completely natural, painless and non-invasive form of therapy intended for all age groups. Differences compared to the devices used in the field of magnetic resonance:

MagnetIQ THERAPY AMPF by selecting the program sets the frequency, the magnetic field density and signal formats. This means that these allowing functions can be targeted (by doctor's orders) effectively used to treat the specific type of pathology.

MagnetIQ THERAPY AMPF excess EM energy is not released in the form of heat (after the currently known therapy and cooling cells causes microdamage to the cell membrane), but in the form of light energy (photons)



Application



PAIN CLINIC

Eliminating or significantly reducing the pain of all etiologies



RHEUMATOLOGY

Rheumatoid arthritis, lumbosacral syndrome, peri-arthritis of shoulder joint, tennis elbow, sudek syndrome



NEUROLOGY

Resolving insomnia and sleep problems



ORTHOPEDY

Bone fractures, sports injuries, chronic pain (bone spurs, ...)



PHYSICAL MEDICINE

Postoperative rehabilitation trauma



SPORTS MEDICINE

All forms of mechanical injuries (muscles, tendons, bones, joints)

Advantages

- A completely natural way of therapy
- Respects the physiological processes in the body
- Painless therapy, not invasive or aggressive Therapy for all age groups
- The treatment is performed in a lying or sitting position
- One treatment lasts 15 minutes
- One therapy is sufficient daily or every other day

Contraindications

- Pacemaker
- Hyperactivity of the thyroid
- Epilepsy
- Extensive internal bleeding
- Shrapnel in the body
- Expectant mothers

Be in charge.

Use Eatelligence.

OrganIQ Food is Eatelligence. It represents a unique plant-based blend made of 100% natural ingredients.

OrganIQ Food has one main goal - to restore the imbalance of the human organism in just three week of daily usage.

Through the period of 21 days, you will feel the general improvement of your body and mind.

To find out how you can get your 3-week system, please visit www.hsttsolutions.com.